BIO SPONGE

COMPRESSAS NEURO-CIRÚRGICAS | NEURO-SURGICAL COMPRESSES

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ENGLISH

Product Description

Biosponge is a neurosurgical pad used in neurological surgery for the purpose of absorbing and removing bodily fluids from the patient and protecting delicate tissues.

Manufacturing Materials

All models of the Biosponge neurosurgical pads are made of rayon and polyester containing barium impregnated strip for visualization of the compress on X-rays. They also contain cotton threads stuck in each pad to remove them during surgery.

Table 1 - Commercial Models

Illustration

<i>Code</i>	<i>Description</i>	<i>Code</i>	<i>Description</i>
2030-08-008	Biosponge 8x8 mm	2030-20-125	Biosponge 20x12,5 mm
2030-10-010	Biosponge 10x10 mm	2030-20-020	Biosponge 20x20 mm
2030-10-025	Biosponge 10x25 mm	2030-20-025	Biosponge 20x25 mm
2030-10-040	Biosponge 10x40 mm	2030-20-040	Biosponge 20x40 mm
2030-12-012	Biosponge 12x12 mm	2030-20-050	Biosponge 20x50 mm
2030-12-025	Biosponge 12x25 mm	2030-20-060	Biosponge 20x60 mm
2030-12-040	Biosponge 12x40 mm	2030-20-070	Biosponge 20x70 mm
2030-12-075	Biosponge 12x75 mm	2030-20-100	Biosponge 20x100 mm
2030-12-512	Biosponge 12,5x12,5 mm	2030-25-025	Biosponge 25x25 mm
2030-12-525	Biosponge 12,5x25 mm	2030-25-050	Biosponge 25x50 mm
2030-12-540	Biosponge 12,5x40 mm	2030-25-070	Biosponge 25x70mm
2030-12-550	Biosponge 12,5x50 mm	2030-25-075	Biosponge 25x75mm
2030-12-570	Biosponge 12,5x70 mm	2030-25-076	Biosponge 25x76mm
2030-12-575	Biosponge 12,5x75 mm	2030-25-100	Biosponge 25x100mm
2030-12-712	Biosponge 12,7x12,7 mm	2030-25-425	Biosponge 25,4x25,4mm

2030-12-725	Biosponge 12,7x25,4 mm	2030-25-476	Biosponge 25,4x76,2mm
2030-12-738	Biosponge 12,7x38,1 mm	2030-28-075	Biosponge 28x75mm
2030-12-776	Biosponge 12,7x76,2 mm	2030-30-030	Biosponge 30x30mm
2030-13-013	Biosponge 13x13 mm	2030-30-070	Biosponge 30x70mm
2030-13-025	Biosponge 13x25 mm	2030-30-090	Biosponge 30x90 mm
2030-13-038	Biosponge 13x38 mm	2030-38-038	Biosponge 38x38 mm
2030-13-051	Biosponge 13x51 mm	2030-38-138	Biosponge 38,1x38,1 mm
2030-13-052	Biosponge 13x52 mm	2030-40-040	Biosponge 40x40 mm
2030-13-060	Biosponge 13x60 mm	2030-40-050	Biosponge 40x50 mm
2030-13-070	Biosponge 13x70 mm	2030-40-070	Biosponge 40x70 mm
2030-13-075	Biosponge 13x75 mm	2030-50-070	Biosponge 50x70 mm
2030-13-076	Biosponge 13x76 mm	2030-50-080	Biosponge 50x80 mm
2030-15-015	Biosponge 15x15 mm	2030-51-075	Biosponge 51x75 mm
2030-15-060	Biosponge 15x60 mm	2030-60-070	Biosponge 60x70 mm
2030-19-019	Biosponge 19x19 mm	2030-70-070	Biosponge 70x70 mm
2030-19-519	Biosponge 19,05x19,05 mm	2030-76-075	Biosponge 76x75 mm
2030-20-010	Biosponge 20x10 mm		

Symbol of the Traceability Label (Label) according to NBRISO15223

Date of Manufacture	Shelf life	Lot Number	Manufacturer Data	Code	Ethylene Oxide Gas
 XXXX	 XXXX-XX	 XXXXX		 XXXX-XX	

TECHNICAL INFORMATION

Indications, Precautions, Restrictions / Warnings, Adverse Effects and Contraindications: Biosponge is indicated for absorption of body fluids in surgical fields, protection of delicate tissues and nerve structures during vacuum suction of fluids released in neurosurgeries. The neurosurgical pad has as its operating principle to absorb and withdraw bodily fluids in neurological surgeries. It must be handled in an appropriate manner and with care. Only trained personnel should handle these products. The surgical pad should be used to absorb blood and other fluids during neurosurgery, protect tissues during vacuum aspiration, and may be used moistened in saline or water. The surgeon must be familiar with and

have sufficient knowledge of the use of the component, surgical technique adopted, precautions and potential risks. Preoperative and surgical procedures, including knowledge of the surgical technique, are important considerations for successful use. Do not use the product if it is damaged. If the product is damaged, in poor condition or with damaged packaging, it should not be used and discarded at the same time. Before closing the surgical field, make sure all neurosurgical pads have been removed. For safety, use X-ray as an aid in this procedure. This product has specific purpose: neurosurgery. Therefore, it is contraindicated for any other use.

Sterilization: This product is supplied

in sterile form, in packaging composed of double blister, Tyvek, Cardboard box and safety seal. Sterilization is done by Ethylene Oxide Gas. Each pack contains 10 units of the same size. If the product has the sterilization validity expired or a damaged packaging, it must be returned to the responsible distributor or directly to Biomechanica.

Product storage and transport: The site must be clean, dry and illuminated in order to maintain the ideal storage and transport conditions, as well as its physical and chemical integrity. The product should be stored and transported in a dry place, with temperature below 30 ° C, relative humidity around 30% to 70% and protected from direct sunlight. It can not be stored directly on the floor (minimum height = 20cm). It can not be placed in very high places, near light bulbs or lamps, which could cause dryness of the packaging or damage to the label. It should not be stored in places where contaminants such as cleaning materials, insectici-

des, pesticides, etc. are stored. Carry carefully to avoid falling and rubbing in order not to damage the packaging. The product should be handled with care to avoid damage that could impair the quality of the material and the safety of the patient. The carrier must be informed about the content and deadline of delivery. The opening of the sealed blister pack can be made manually, in a suitable place (surgical environment), by pulling the Tyvek. This opening should be performed only in the surgical act to avoid contamination of the product and by specialized persons. Carry carefully to avoid falling and friction in order not to damage the package and the component, thereby ensuring the sterility of the material.

Identification: Biosponge – Neurosurgical Pad has labels that include the information: Manufacturer Data; Technical manager; Registration Number in ANVISA; Product Technical Name; Product Commercial Name; Commercial Product Model; Raw Material; Number

of items in the package; Lot Number; Product Code; Manufacture Date; Expiring Date; Manufacturer's Address.

Customer Service: If there is a need to make a complaint regarding the use of Biosponge related to any adverse effect that affects the safety of the user, the surgeon responsible should communicate this adverse event to the competent sanitary body and Biomechanica through the email tecnovigilancia@biomecanica.com.br or by telephone +55 14 2104 7900. In case of doubts the surgeon in charge or the health professional can communicate the adverse event through the System of Notifications in Sanitary Surveillance in the ANVISA website: <http://www.anvisa.gov.br/hotsite/notivisa/index.htm>

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